

INSTALLATION GUIDE SUCTION SUPPORT BAR

IMPORTANT

- This support bar is designed to assist with balance, and not suitable to support the entire body weight of the user.
- This support bar must not be installed on surfaces that are: dirty, greasy or wet; porous (e.g. unfinished wood); uneven or not flat.
- This support bar must not be installed over grout lines or other gaps in the wall.
- Check that the support bar is fully secure prior to each use.

GLOW-IN-THE-DARK MODELS

To achieve sufficient luminance, the support bar must be exposed to ultra violet light, such as the sun or a fluorescent light source.

CARE AND MAINTENANCE

- Clean with a soft cloth and warm soapy water.
- Do not use harsh detergents, corrosive chemicals, or abrasive cleaners.



1. Clean and dry the wall and suction pads.



2. Lift up (unlock) the blue tabs on the support bar.




3. Position the support bar into the desired position, ensuring the pads are not sitting on top of any grout lines or gaps in the wall.



4. Push the support bar firmly against the wall and hold this position while performing the next step.



5. Push the tabs the whole way down.



6. Check that the support bar is fully secure before each use.